

Finished Product Inspection Report

Report NO: COVID-19Ag-cp-ST210403

Product name	COVID-19 (SARS-CoV-2) Antigen test kit (kolloidales Gold)		Specification	Cassette, 1 pcs/foil pouch
LOT	ST210403		Quantities	2000000 pcs
Applying date	2021/04/23		Test date	2021/04/23
Report time	2021/04/23			
NO	Test Item	Standard code		Test results
1	Physical character	Appearance	The test should be clean and complete, no burr, no damage and non-pollution. The shell of the test cassette should be flat, the upper and lower covers should be evenly closed, and there should be no obvious gap. The inner test strip should be firmly attached without wobble. The diluent should be free of foreign matter.	Pass
		Strip width	Should be not less than 2.5mm	3.0mm
		Rate of migration of liquid	Should be not less than 10mm/min	40.3-45.9 mm/min
2	Minimum detection limit	The enterprise reference products S1, S2 and S3 shall be tested, wherein S1 shall be negative and S2 and S3 shall be positive.		Pass
4	Positive compliance rate	5 copies of positive reference products of the test company shall be all positive, with a positive compliance rate of 100%.		Pass
5	Negative compliance rate	5 copies of negative reference products of the test company shall be all negative, with a negative compliance rate of 100%.		Pass
6	Repeatability	Test 1 copy of enterprise positive reference, test at least 10 times, are positive.		Pass
Inspection standard	According to The enterprise standard SLSCQJ-ZL-163			
Inspection conclusion	This product has been tested, it is suitable for The enterprise standard SLSCQJ-ZL-163			

Quality director: *Jing Xu*

Checker: *Xinyue peng*



Inspector: *Li Zhong*